


Home Phototherapy

Model: 027



## Technical Specification

- 4 \* 20 watt Philips lamps (Code no. TL 20W/52)
- Timer for lamp life time measuring
- Timer for treatment duration measuring
- Cooling system for electronics and electric circuit
- 4 folding leg

## Alarms

- Audible alarms by terminality treatment time
- Visible alarm by terminality lamp life time

## Electrical requirement

- 220V AC, 50 Hz, 100W

## Dimension and weight

- Equipment :
  - Length :68
  - Width :50
  - Height: 9
  - Weight: 7.5 Kg
- Bag:
  - Length: 71
  - Width: 58
  - Height: 13
  - Weight: 5Kg