


Infant Incubator

توسان Tosan

Manufacture & Supply of Neonatal System



Model: 111



Technical Specification

Display:

Digital air temperature thermometer

- Resolution : 0.1° C
- Range : 28 to 40° C
- Controllable range : 28 to 37° C override to 39° C

Digital Skin temperature thermometer

- Resolution : 0.1° C
- Range : 33 to 38° C
- Controllable range : 33 to 37° C

Heater - power indicator

- Flashing L.E.D

Performance

- Internal noise level : less than 48 dB
- Warm-up time : up to 10° C 30 minutes
- Power consumption : 350 W
- Humidification : up to 90 % RH
- Controlable range : 30 - 90
- Disposable filter

Alarms

- Air sensor fault
- Skin sensor fault
- Air circulation failure
- Power failure
- Air set deviation +/- 2° C
- Low battery
- Low water

Electrical requirements

- Power supply :220V 50Hz 3A

Dimension and weight

- Mattress : 34 * 60 cm
- Incubator length , width , height :92,63,145
(Adjustable from 125 to 145 cm if
equipped with electrical jack)
- Weight : 90 Kg

Available options

- Iv pole
- Cabinet
- Height adjustment

Head Office :No,19.Shabnam (west sarv)St.

Sarve Square , Saadatabbad , Tehran , IRAN

Tel:(+98-21)22072192

Fax:(+98-21)22362814

www.tosanco.com

Email:info@tosanco.com