



Radiant warmer model 052

CE
1254

Warmer is used to locally warm up the neonate's body

Technical Specification:

Digital skin temperature thermometer:

- Skin temperature display
- Resolution: 0, 1 °C
- Range: 28 to 45 °C
- Controllable range: 28 to 37 °C (manual and skin mode)
- One 20 Watt halogen lamp.
- Adjustability of the elevation
- With 4 wheels to facilitate to move the heater
- Warning alarm to notify the nurse every 10 minutes.

Technical specifications:

- Power source: 220 volts, 50 Hz, 3Amper
- Dimensions: length: 52, width: 78, height: 163 cm



Head office; No.14 .Shabnam (west sarv) St.
Sarve Square, Saadatabbad ,Tehran, Iran.
Tel: +98 21 2072192
Fax: +98 21 2362814
www.tosanco.com
Email:info@tosanco.com

Tosan co. Continuously improves the productions and for this reason the specifications of the purchased equipments may be optimized comparing to those explained in this brochure.