

# Neonatal Resuscitator

توسان Tosan

Manufacture & Supply of Neonatal System



Model: 083



1254



## Technical Specification

### Display:

- Skin temperature thermometer
- Resolution:0.1°C
- Temperature range : 28° to 39.3°C
- Controllable temperature range : 28° to 37°C  
(Manual and Skin mode)

### Heater power indicator

- 10% to 100% LED Power

### Performance

- Adjustable timer
- 2 \* 35 Watt halogen Lamps
- Mechanical tray tilting up 12 degree front and rear
- Electromechanical tray height adjustment  
from 80 to 100 cm
- IV stand
- Monitor tray (option)
- Radiology cassette
- 2 shelf
- Rotary Head while radiology
- Adjustable Oxygen flow
- Adjustable Suction vacuum

### Alarms

- Skin sensor fault
- Skin temperature overheat
- Temperature set deviation
- Power failure

### Electrical requirement

- Power supply : 220 V 50 Hz 3A

### Dimension and weight

- Mattress : 34 \* 60 cm
- Length,width,Height :99,66,198
- Weight :