



Phototherapy unit model 023

Technical Specifications:

- 4 x 20 watt Philips lamps (Code No. TL 20W/52).
- Timer for lamp life time measuring
- Timer for treatment duration measuring
- Height adjustable
- Adjustable angle of light exposure.
- Usable over the incubator or cot

Electrical Requirements:

- Power supply: 220V, 50Hz, 3A.

Dimensions & weight:

- Length: 70cm
- Width: 66cm
- Height: 165 cm
- Weight: 21.3 kg

Head office; No.14 .Shabnam (west sarv) St. Sarv Square, Saadatabbad ,Tehran, Iran.

Tel: +98 21 2072192 Fax: +98 21 2362814

www.tosanco.com Email:info@tosanco.com